medipal

ALCOHOL
DISINFECTANT
Wipes

With
TRITEX®
TECHNOLOGY

200
WIPES

For disinfecting surfaces of non invasive medical devices

Medipal Alcohol Disinfectant Wipes: are intended for disinfecting surfaces of non invasive medical devices and equipment within healthcare environments. **Professional Use Only.** Refer to training documentation before use.

DIRECTIONS FOR USE: Use PPE as required. Use one wipe to disinfect, holding the wipe flat and wiping the surface in an S-shape ensuring the entire surface is wet. Use a fresh wipe if your wipe becomes soiled or dry. For visibly soiled surfaces, use one wipe to remove visible soil prior to disinfection. Do not use on surfaces sensitive to alcohol. Allow the surface to dry naturally in line with stated contact time. Do not re-use wipes. Do not use on patients or invasive medical devices. This product is not intended for end point sterilisation of medical devices.

STORAGE: Keep pack tightly closed when not in use. Store in a cool, dry & well ventilated place. Once opened use within one month.

DISPOSAL: Dispose of wipes and packaging in accordance to local legislation. Do not flush wipes down the toilet. Do not macerate.

For further information about specific microorganisms please refer to the technical data sheet.

Microorganism	Contact Time
Tuberculosis/Mycobacteria	
EN14348	1 min
Bacteria	
EN13727	1 min
EN16615	1 min
Fungi/Yeast	
EN13624 (yeastcidal)	3 min
EN16615 (yeastcidal)	3 min
Virus	
EN14476 (Enveloped viruses)	1 min

Available from
www.palinternational.com

**Medipal Alcohol Disinfectant Wipes
Danger**

Ingredients/content at manufacture:
Isopropanol (CAS:67-63-0) 70% v/v (629g/kg)
H225: Highly flammable liquid and vapour. H319: Causes serious eye irritation. H336: May cause drowsiness or dizziness. EUH066: Repeated exposure may cause skin dryness or cracking.
Safety: P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P261: Avoid breathing vapours. P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing. P337+313: If eye irritation persists get medical advice/attention. P233: Keep container tightly closed. P403 + P235: Store in a well-ventilated place. Keep cool.

In accordance with (EC) No 1272/2008
Keep out of reach and sight of children.

SAFETY:

STORAGE:

DO NOT
STORE IN
DIRECT
SUNLIGHT

EU Importer:

Pal Hygiene Products Limited,
Unit 5B & Unit 5H, Fingal Bay
Business Park, Balbriggan, Co. Dublin.
Republic of Ireland, K32 NY57
T: +353 15060172
E: importer@palinternational.com

Manufactured in the UK by:

Pal International Limited
Unit 3 Mountpark, Bardon II,
Victoria Road, Ellistown,
Coalville, LE67 1FA
T: +44 (0) 1455 555700
E: info@palinternational.com
W: www.palinternational.com

Advena Limited
Tower Business Centre
2nd Floor.
Tower Street
Swatar, BKR 4013
Malta
E: info@advena.mt

For batch
number
and use-by
date see
side of
package

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