ANIOSYME5 SYNERGY

5 ENZYMES FOR 1 UNIQUE PRODUCT

LABORATOIRES

ANIOS

An Ecolab Company

5L

5C2235861005LEMR-230330
3131162F/30MAR23/50

ANIOSYME5

SYNERGY

UFI: 6UGQ-8Q3E-YF0R-YAEG

**ATTENTION • WARNING • ACHTUNG • WAAR-
SCHUWING • ATTENZIONE • ATENÇÃO • ATENCIÓN •
ΠΡΟΣΟΧΗ • ΔΙΚΤΑΤ**

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EN 5 enzymes cleaning product for instrumentation

INDICATIONS: Cleaning of surgical and medical equipment, endoscopes and medical devices, in soaking, ultrasonic tanks, automatic washing machine and washing tunnel. **COMPOSITION**: Non-ionic and cationic surfactants, chelating agent, stabilizing agent, enzymatic complex, excipients. **INSTRUCTIONS FOR USE**: Prime the associated pump before first use. **Manual use**: Concentrated solution. Use at a dilution from 0.5% to 2%. Effective from 5 minutes. Adapt the concentration and the contact time of the dilution water to the degree of soiling. Brush if necessary. Rinse thoroughly. Maximum temperature: 35°C. **Automated use**: Concentrated solution. Use at a dilution from 0.1% to 0.5%. Effective from 5 minutes. Adapt the concentration, contact time and temperature of the dilution water to the degree of soiling. Rinse thoroughly. Maximum temperature: 60°C. **PRECAUTIONS FOR USE**: (P273) Avoid release to the environment. (P280) Wear eye protection/face protection. Dispose of the material and its container by following all the warnings of use. Reserved for professional use. Medical Device. Made in France.

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