ANIOSPRAY QUICK

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 LABORATOIRES
ANIOS
An Ecolab Company

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EN INDICATIONS: Fast disinfection of previously cleaned, non-immersible, non-invasive, alcohol-resistant surfaces, medical equipment and medical devices; and general surfaces (biocide application). **COMPOSITION:** 100 g contains : 55 g Ethanol (CAS N° 64-17-5), 0.11 g DMPAP, denaturing agent, perfume. **MICROBIOLOGICAL PROPERTIES:** Bactericidal in 1 min: EN 13727* (30 sec), EN 17387* (1 min). Yeastcidal according to EN 13624*/EN 17387* in 30 sec. Fungicidal according to EN 13624*/EN 17387* in 5 min. Tuberculocidal according to EN 14348* in 30 sec. Mycobactericidal according to EN 14348* in 5 min. Active against enveloped viruses according to EN 14476*/EN 16777* in 30 sec. Limited spectrum virucidal according to EN 14476*/EN 16777* in 5 min (Adenovirus in 30 sec / Norovirus (MNV) in 5 min). * = clean conditions. **INSTRUCTIONS FOR USE:** Ready-to-use solution. Use on previously cleaned surfaces and medical devices. Spray on the areas to be treated until they are totally moistened. Apply in sufficient quantity (+/- 40 ml/m²). Complete wetting of all surfaces to be disinfected is important. Respect the indicated contact time for the required antimicrobial activity. Rinsing is not necessary unless the surface comes into contact with the skin or mucosa. For the frequency of use, refer to the hygiene plan in place. Please refer to the manufacturer's instructions for use for cleaning and disinfecting medical devices. Remarks: Before application, check the compatibility of the areas to be treated with alcohol. Electrical devices being disinfected must be switched off. Scientific dossier (efficacy, compatibility, etc.) available on request. **PRECAUTIONS FOR USE:** (H225) Highly flammable liquid and vapour. (H319) Causes serious eye irritation. (P210) Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. (P280) Wear eye/face protective equipment. Dispose of the product and its container, following all precautions for use. Biocidal product intended for the disinfection of surfaces (Group 1 - PT2). Reserved for professional use. (a) Keep out of the reach of children. Anti-Poison Center Tel: 210 7793777 (Greece), Registration number: 47438/24. (cy) 1401 (Cyprus). **EL ΕΝΔΕΙΞΕΙΣ:** Ταχείας δράσης απολύμανση επιφανειών, ιατρικού υλικού, καθώς και μη βυθιζόμενων και μη επεμβατικών ιατρικών συσκευών,

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