ANIOSYME XL3

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Centre antipoison ☎ +32(0)70-245245 ☎ +352(0)8002-5500.

EN Detergent Pre-disinfectant for instrumentation

INDICATIONS: Cleaning and disinfection of medical devices: surgical and medical instruments, heat-sensitive instruments and endoscopy equipment, before sterilisation.

COMPOSITION: Quaternary ammonium carbonate, non-ionic surfactants, enzyme complex (protease, amylase and mannanase), chelating agent, perfume, colouring, excipients.

MICROBIOLOGICAL PROPERTIES: Bactericidal* according to EN 13727/EN 14561 in 5 min, 0.5%. Yeastocidal* according to EN 13624/EN 14562 in 5 min, 0.5%. Tuberculocidal* according to EN 14348/EN 14563 in 30 min, 1%. Active on enveloped viruses* according to EN 14476/EN 17111 in 15 min, 0.5%. Disinfection of instruments: according to VAH (bactericidal and yeastocidal)* in 15 min, 0.5%. * = dirty conditions. **INSTRUCTIONS FOR USE:**

Concentrated solution. Prime the associated pump before first use. To be used in soaking at the required dilution, respecting the contact time. Totally immerse the medical device. Keep the soaking bath covered during the contact time. Brush if visible soiling is present. Thoroughly rinse treated medical devices with tap water to remove all traces of product and organic soiling. Replace the solution as soon as it is visibly soiled and at least once per day to maintain the effectiveness of the product. Scientific dossier (efficacy, compatibility, etc.) available on request. **PRECAUTIONS FOR USE:** (P273) Avoid release to the environment. (P280) Wear protective gloves/ eye protection/ face protection. (P305 + P351 + P338) IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. (P310) Immediately call a POISON CENTER or doctor/physician. Dispose of the product and its container, following all precautions for use. Reserved for professional use.

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