SURFA'SAFE *Premium*

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GR24190750MCAR-210823

SURFA'SAFE *Premium*

EN Disinfecting detergent solution for surfaces and

medical devices. INDICATIONS: Disinfectant detergent solution for cleaning and disinfecting non-invasive, non-immersible medical devices. Biocide application: for general surfaces.

COMPOSITION: 100 g contains: 0,3 g of Didecyltrimethylammonium chloride (DDAC) (CAS No. 7173-51-5), non-ionic surfactants, chelating agent, excipients. **MICROBIOLOGICAL PROPERTIES:** Bactericidal in 2 min: EN 13727* (1 min), EN 16615* (2 min). Yeastcidal in 2 min: EN 13624* (2 min), EN 16615* (2 min). Tuberculocidal in 60 min: EN 14348*. Active against enveloped viruses according to EN 14476* in 1 min. Active in accordance with EN 14476* against Polyomavirus SV 40 in 5 min. Active in accordance with EN 14476* against Rotavirus in 30 sec. Limited spectrum virucidal according to EN 14476* in 30 min (Adenovirus in 30 min / Norovirus (MNV) in 15 min). *=dirty conditions.

INSTRUCTIONS FOR USE: Ready-to-use solution. Spray the disinfectant detergent solution on a non-woven cloth and wipe. Spread evenly the product (+/- 6 ml/m² or around 10 sprays). Respect the indicated contact time for the required antimicrobial activity. Rinsing is unnecessary except for the areas intended to come into contact with the skin or mucosa. Please refer to the manufacturer's instructions for use for cleaning and disinfecting medical devices. For the frequency of use and cleaning of the application equipment, refer to the hygiene plan in place. Scientific dossier (efficacy, compatibility, etc.) available on request. **PRECAUTIONS FOR USE:** Dispose of the product and its container, following all precautions for use. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State. Biocidal product intended for the disinfection of surfaces (Group 1 - PT2). Medical Device. Reserved for professional use. Made in France.

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SURFA'SAFE Premium

EN Disinfecting detergent solution for surfaces and medical devices. **INDICATIONS:** Disinfectant detergent solution for cleaning and disinfecting non-invasive, non-immersible medical devices. Biocide application: for general surfaces.

COMPOSITION: 100 g contains: 0,3 g of Didecylidimethylammonium chloride (DDAC) (CAS No. 7173-51-5), non-ionic surfactants, chelating agent, excipients. **MICROBIOLOGICAL PROPERTIES:** Bactericidal in 2 min: EN 13727* (1 min), EN 16615* (2 min). Yeastcidal in 2 min: EN 13624* (2 min), EN 16615* (2 min). Tuberculocidal in 60 min: EN 14348*.

Active against enveloped viruses according to EN 14476* in 1 min. Active in accordance with EN 14476* against Polyomavirus SV 40 in 5 min. Active in accordance with EN 14476* against Rotavirus in 30 sec. Limited spectrum virucidal according to EN 14476* in 30 min (Adenovirus in 30 min / Norovirus (MNV) in 15 min). *dirty conditions.

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+43 1 715 25 50-0 / ECOLAB GmbH -