SURFANIOS *Premium*

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EN INDICATIONS: Cleaning and disinfection of non-invasive medical devices prior to use and according to their instructions for use. **COMPOSITION:** N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine (CAS N° 2372-82-9: 51 mg/g), Didecyldimethylammonium chloride (CAS N° 7173-51-5: 25 mg/g), non-ionic surfactants, chelates, perfume, colorant and excipients. **MICROBIOLOGICAL PROPERTIES:** Bactericidal in 15 min, 0,25% : EN 13727* (5 min, 0,25%), EN 16615* (15 min, 0,25%). Yeastcidal in 15 min, 0,25% : EN 13624* (15 min, 0,25%), EN 16615* (15 min, 0,25%). Tuberculocidal in 30 min, 0,5% : EN 14348*. Active against enveloped viruses according to EN 14476* in 5 min, 0,25%. Limited spectrum virucidal according to EN 14476* in 30 min, 2% (Adenovirus in 30 min, 1% / Norovirus (MNV) in 30 min, 2% or 60 min, 1%). *dirty controls. **INSTRUCTIONS FOR USE:** Concentrated solution. Dilute the concentrate with water to the desired concentration. To dilute e.g. to a concentration of 0,25% (2,5 ml/l), add 20 ml for 8 litres of cold or hot water (max. + 60 °C). Apply in sufficient quantity (+/- 30 ml/m²). Respect the indicated contact time for the required antimicrobial activity. Rinse thoroughly with tap water to remove any trace of the product. For the frequency of use and cleaning of the application equipment, refer to the hygiene plan in place. Scientific dossier (efficacy, compatibility, etc.) available on request.

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EN PRECAUTIONS FOR USE: Contains N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine and didecyldimethylammonium chloride. (H290) May be corrosive to metals. (H315) Causes skin irritation. (H318) Causes serious eye damage. (H410) Very toxic to aquatic life with long lasting effects. (P273) Avoid release to the environment. (P280) Wear protective gloves/ eye protection/ face protection. (P305 + P338) IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. (P310) Immediately call a POISON CENTER or doctor/physician. Dispose of the material and its container by following all the warnings of use. Keep out of the reach of children. Anti-Poison Center Tel: 210 7793777 (Greece), 1401 (Cyprus). Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State. Reserved for professional use. Made in France.

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